A photograph of a male chiropractor in a light blue polo shirt and glasses, wearing a white glove on his right hand, examining a female patient's neck. The patient is lying on a red chiropractic table with her eyes closed. The background is a plain wall.

Ultimate  
Reliability

MANUFACTURED & MARKETING BY: 

**ORTHOCURE HEALTHCARE PRIVATE LIMITED,**

Plot. No. 188, J-2, Udyog Vihar,

Phase IV, Sec 18, Gurugram, Haryana, 122015

Customer Support on : [support@orthocure.co.in](mailto:support@orthocure.co.in)

For Purchase Contact : 7262001001 : [sales@orthocure.co.in](mailto:sales@orthocure.co.in)

**CHIROPRACTIC**

**PHYSIOTHERAPY**

**POSTURE & FITNESS PRODUCTS**

Website : [www.orthocure.co.in](http://www.orthocure.co.in)



Products | 2022  
Catalogue | 2023

# Back Support



## SPECIFICATIONS


- Colour : Black
- Package : 37.69 x16.1 x9.5 cm
- Weight : 0.181.44 gm

## FEATURES

- It is recommended by chiropractors, orthopaedist & physiotherapists to prevent neck, back, hip, leg and sciatic nerve pain with regular use
- It Improves blood circulation
- Reduce muscle Fatigue and tension
- Improves your working efficiency
- Maintains the natural curvature of the spine



# Cervical Pillow



## SPECIFICATIONS

- Weight : 1.225 kg
- High quality memory foam
- Hollow head resting area
- Neck support area
- Ergonomic design

## FEATURES

- Orthocure medium density memory foam pillow provides maximum Comfort
- Use this pillow for regular use below your neck/head while sleeping, resting or use it as directed by your physiotherapist/physician
- This pillow is more comfortable then ordinary pillow to sleep with



# Adjuster



## SPECIFICATIONS

- Dimension : 7.6 x 3.5 x 0.75 Inch
- Weight : 260 gm
- Colour : Blue & Red
- Adjustable Force : 50-100 N
- Package Included : Single stylus, T-shape stylus, U-Shape stylus


## FEATURES

- Activator is a spring-loaded hand-held instrument to provide quick and low-force impulse to specific areas
- It improve joint pain and spasm : Nack pain, Back pain
- Use full in adjusting stiff vertedration syndrome, movement disorders, thoracic spine syndrome

Launching on 1st Feb 2022



# Peripheral / Extremity Drop



## SPECIFICATIONS

- Wooden Base with a firm cushion
- 8" x 6" x 2.5" (20x15x6 cm)
- Weight : 800gm (Approx.)
- Colour : Black
- Drop Height : 8.4mm (Maximizes velocity)
- With Strap Handle

## FEATURES

- The peripheral drop has a flat surface with straight drop
- It is ideal for adjusting peripheral extremities
- This can help in improving mobility of the affected joints
- Helps in manipulation techniques which can not be done manually
- Safe to use on patients & easy to carry



# Lumbo-Pelvic Drop



## SPECIFICATIONS


- Size : 51 x 31 x 9 cm
- Weight : 4 kg
- Drop height : 14 mm
- Tapers to height of 4.4 cm at hinged end
- Colour : Black

## FEATURES

- Lumbo pelvic drop is ideal for Pelvic, Lumbar, thoracic & extremity adjusting
- Herniated discs, pain in the hips and lower back, pelvic girdle dysfunction
- With a quick & crisp drop, It tapers toward the inferior end for patient comfort



# Cervical Drop



## SPECIFICATIONS

- Headpiece : 33 x 25 x 10 cm
- Durable base : 30 x 33 cm
- Weight : 3 kg
- Drop Height : 16 mm

## FEATURES

- With both straight & forward motion drops, you can utilize it for a variety of cervical adjusting techniques including side posture toggle recoil
- The stabilization base distributes the weight and force of cervical adjusting over a broad area
- The contoured face cushions is provided for patients comfort



## Pevic Wedges



### SPECIFICATIONS

- Approx 20 x 12 x 11 cm
- Weight : 0.7 kg (Approx)
- Wooden Base with a firm cushion
- Colour : black
- Solid wood interior
- High density from

### FEATURES

- Bottom of each block is covered in vinyl & secured with upholstery nails
- These wedges are placed under each side of the pelvis thus drawing your disc away from the nerve being affected
- Safe to use on patients & easy to carry

## SI Belt



### SPECIFICATIONS

- Available in four different sizes : S, M, L, XL
- Made with breathable and hypoallergenic material for comfort and durability
- Superior compression Elastic helps reduce muscle spasms and encourage correct posture

### FEATURES

- It can be used in various conditions like back pain, hip joint pain, sciatica, SI joint, hip degeneration, disc herniation, facet syndrome
- It can be used in pubic symphysis dysfunction pelvic girdle dysfunction & piriformis syndrome
- It is a 3 inch wide belt which relieves the stress on the joint and provides stability

## Rehab Wall Station



### SPECIFICATIONS

- 200 x 22. 61 x 15. 24 cm
- Weight : 12 kg
- High quality anodized Aluminum construction
- Protatable Inclinometer

### FEATURES

- The Wall Station is a compact total body rehabilitation system designed for in clinic strength training.
- It features the familiar colors & resistance levels of clip connect tubing.
- Upper & lower extremities features a multidimensional design with three planes of moment.
- Inclinometer helps to check passive range of motion.